

INTRODUCTION

Prostate-specific membrane antigen (PSMA) PET/CT imaging has provided significant advances in diagnosing and treating patients with prostate cancer. PSMA PET/CT has been used in AI models to evaluate lesions; however, few studies have assessed the transferability of AI models across the commonly used tracer types.

In this analysis, we deployed a CNN-based AI model designed to identify lesions on both ⁶⁸Ga-PSMA-11 PET/CT and ¹⁸F-DCFPyL PET/CT images and evaluated the generalizability and performance on ¹⁸F-PSMA-1007 PET/CT images, which has a different biodistribution.

MATERIALS AND METHODS

- A set of N=169 images, comprising ⁶⁸Ga-PSMA-11 PET/CT images from 89 patients (126 images, 1-2 imaging time-points/patient) and ¹⁸F-DCFPyL PET/CT images from 43 patients (43 images, 1 imaging time-point/patient), all diagnosed with metastatic prostate cancer, were used for training a lesion detection CNN employing a retina U-net architecture.




Image from Jaeger et al., 2018

Figure 1: Architecture for the Retina U-net [1]

- The CNN's performance was assessed separately in external validation datasets with ⁶⁸Ga-PSMA-11 PET/CT images (200 patients, 337 images), ¹⁸F-DCFPyL PET/CT images (27 patients, 27 images) and an additional dataset of ¹⁸F-PSMA-1007 PET/CT images (13 patients, 18 images).
- Lesions were manually delineated for comparison and reviewed by a clinician with extensive training in reading PSMA PET/CT.
- Detection performance was evaluated using the sensitivity and the number of false positives per image (FPs/image).
- Performance was summarized for all lesions, for lesions with $SUV_{max} > 2.5$ g/ml, and for lesions with $SUV_{max} > 4$ g/ml.

RESULTS

- In the ⁶⁸Ga-PSMA-11 test images (N= 337), an average of 4 lesions per image (range: 1-100 lesions per image) were delineated. The median detection sensitivity for all lesions was 92% (IQR: 43%-100%) with 1 (IQR: 0-2) FPs/image




Figure 2: Comparison of the Lesion Detection Sensitivity (%) and False Positives (N) across different tracers

- For the ¹⁸F-DCFPyL test images (N= 27), an average of 20 lesions per image (range: 1-157 lesions per image) were delineated. The median detection sensitivity for all lesions was 100% (IQR: 82%-100%) with 2 (IQR: 1-2.5) FPs/image.
- Similarly, for the ¹⁸F-PSMA-1007 test images (N=18), an average of 7 lesions per image (range: 1-29 lesions per image) were delineated. The median detection sensitivity for all lesions was 100% (IQR: 95%-100%) with 2 (IQR: 1-4) FPs/image.




Figure 3: Patient-level detection performance with two example patients illustrating the median performance and false positive distribution across the different tracers: A) ⁶⁸Ga-PSMA-11, B) ¹⁸F-DCFPyL, C) ¹⁸F-PSMA-1007

KEY FINDINGS

A CNN detection model, trained on ⁶⁸Ga-PSMA-11 and ¹⁸F-DCFPyL PET/CT images, demonstrated high median sensitivity and minimal false positive rates in the ⁶⁸Ga-PSMA-11, and ¹⁸F-DCFPyL external holdout datasets, and on the ¹⁸F-PSMA-1007 images.

Although the ¹⁸F-PSMA-1007 dataset was small, this preliminary data shows the potential for generalizability across different PSMA tracers, including those with different biodistributions such as ¹⁸F-PSMA-1007.

Investigating generalizability in larger, diverse cohorts is vital to fully understand the implications for model stability and clinical reliability. Analysis is ongoing to include more ¹⁸F-PSMA-1007 images in this study.

REFERENCES

- [1] Jaeger, P. F. et al. Retina U-Net: Embarrassingly Simple Exploitation of Segmentation Supervision for Medical Object Detection. *Proc. Mach. Learn. Res. NeurIPS 2019* 1–12 (2018).

DISCLOSURES

Authors OL, BD and TGP are employed by AIQ Solutions. RF is a scientific advisory board member of AIQ Solutions. MD is an AIQ Research Fellow. Remaining authors have no relevant relationships to disclose.

CONTACT INFORMATION

✉ ojaswita.lokre@aiq-solutions.com